




# Red Therapy Light

*Deep Red •660nm-50%+Near Infrared 850nm-50%*

- ✓ Skin repair of your acne on the face, whitening and scar
- ✓ Pain relief for shoulder, waist and abdomen
- ✓ Promote muscle healing wound healing
- ✓ Improve Sleep Blood Circulation



# The Benefits of Red Light Therapy



Low power  
Enhanced  
microvascular




Promote blood  
circulation, increase  
oxygen consumption  
accelerate detoxification



Improve blood  
& lymphatic  
system



Stimulating  
cell




Re-new  
Collagen



Repair cell  
Healthy skin

# Red Light Therapy



**Benefits of Red Wavelength**



### Reduction of lines, wrinkles, and hyperpigmentation:

660nm can support the production of collagen, assisting with plumper, more radiant-looking skin, and a more uniform texture.

### The treatment of acne:

The 660nm wavelength has also been shown to help reduce inflammation and breakouts.



### An improvement in skin disorders such as psoriasis:

One study showed that the use of the 660nm wavelength used in conjunction with hematoporphyrin derivative (HPD) helped to eradicate the symptoms associated with psoriasis within 17 days.

### Reduced training fatigue:

The 660nm wavelength teamed with the 830nm near-infrared wavelength has been proven to delay the development of fatigue in the muscles, and enhance skeletal muscle performance. (Athletes, take note.)




### Reduced neuropathic pain:

A study investigating sciatic nerve pain in rats found that the 660 nm wavelength significantly helped to reduce pain.

### Hair regrowth:

Red light therapy containing the 660nm wavelength has been proven to be a safe and effective method of stimulating hair growth in both men and women. One study hypothesized that the light stimulated stem cells in the hair follicle, encouraging the hair follicle into an anagen or growth phase.



# Infrared Light Therapy



**Benefits of NIR Wavelength**



### Anti-inflammatory benefits:

850nm wavelengths can help to reduce joint and muscle pain and diminish general inflammation in the body.

### Enhanced muscle recovery:

A study observed the use of 850nm wavelengths on athletes, and found that usage of the near-infrared light increased muscle mass after training, and decreased inflammation and oxidative stress in muscle biopsies.



### Healing of wounds in the skin:

Lesions in the skin heal faster when exposed to 850nm wavelength light therapy.

### Improved bone healing:

The 850nm wavelength encourages resorption and formation in the bone cells around the location where repair is needed, without causing any change to the bone structure.




### Orthodontics and tooth alignment:


The use of the 850nm wavelength can foster the speedy realignment of teeth for patients undergoing orthodontics.

# Product Specifications


## AG-300 Plus




LED Power :300W  
 Dimensions: 15.7" x6.7" x3"  
 Power Consumption: 90W ± 5%  
 Amps @110v:0.82A  
 Amps @220v: 0.41A  
 Treatment Area: See Chart  
 LEDs:60 x 5W  
 Focusing Lens Beam Angle :  
 60 Degree  
 Cooling fans:1 Fan  
 Switchable: YES  
 Weight:3.5KG  
 Warranty:2 Years  
 Lifespan:100,000 Hrs  
 EMF Emission:0.0 μT @ 6"  
 Best For:Targeted Treatment  
 Max Area:64" x57"



# AG-480 Plus




LED Power :480W  
 Dimensions: 19" x9" x3"  
 Power Consumption: 140W ± 5%  
 Amps @110v:1.27A  
 Amps @220v: 0.64A  
 Treatment Area: See Chart  
 LEDs:96 x 5W  
 Focusing Lens Beam Angle :  
 60 Degree  
 Cooling fans:2 Fans  
 Timable: YES  
 Weight:4.0KG  
 Warranty:2 Years  
 Lifespan:100,000 Hrs  
 EMF Emission:0.0 μT @ 6"  
 Best For:Targeted Treatment  
 Max Area:64" x58"







## AG-900 Plus




LED Power :900W  
 Dimensions: 36" x9" x3"  
 Power Consumption: 240W ± 5%  
 Amps @110v:2.18A  
 Amps @220v: 1.10A  
 Treatment Area: See Chart  
 LEDs:180x 5W  
 Focusing Lens Beam Angle :  
 60 Degree  
 Cooling fans:3 Fans  
 Timable: YES  
 Weight:7.5KG  
 Warranty:2 Years  
 Lifespan:100,000 Hrs  
 EMF Emission:0.0 μT @ 6"  
 Best For:Targeted Treatment  
 Max Area:76" x54"



## AG-1485 Plus



LED Power :1485W  
 Dimensions: 36" x11.4" x3"  
 Power Consumption: 400W  $\pm$  5%  
 Amps @110v:3.64A  
 Amps @220v:1.82A  
 Treatment Area: See Chart  
 LEDs:297 x 5W  
 Focusing Lens Beam Angle :  
 60 Degree  
 Cooling fans:4 Fans  
 Timmable: Yes  
 Weight:14.0KG  
 Warranty:2 Years  
 Lifespan:100,000 Hrs  
 EMF Emission:0.0  $\mu$ T @ 6"  
 Best For:Targeted Treatment  
 Max Area:76" x60"



### Instructions:

1. *Distance: Apply the light from 0 - 30 cm to the body.*
2. *Session Time: Exposure each area for 10 - 25 minutes, depending on distance.*
3. *Repeat on other areas if desired.*
4. *Frequency: Do this 2 times to 14 times a week until results are felt.*
5. *Maintenance: Continue using several times a week.*
6. *Light Area: Led Light can be applied to any area of the body, except the eyes*


### Note:


Indoor use only : Light is not waterproof , can not used in water or in wet places.  
 Repairing : If your light need after-sale service , please contact us directly .  
 Installation : To avoid the danger of electric shock, please turn off the powerswitch before installation.  
 Extremely bright : Do not look directly at the LED light when it's works.

# More Detail

## Secondary Focusing Lenses

Not only are we able to emit the light in an even manner, we're able to control the angle of this beam of light to focus it tighter with more centralized intensity or wider for more coverage in the exact same manner as a focusing flashlight.





## Red indicator light

Press the “Red” button, the “Red” indicator light turn on, Red wavelength and fans running: Press “up” button, to set time, each change with 5 minutes increase to 60 minutes: Press “down” to set time, each change with 5 minutes decrease to 0..




## Green indicator light

Press the “NIR” button, the “Green” indicator light turn on, N. IR wavelength and fans running; Press “up” button to set time each change with 5 minutes increase to 60 minutes: Press “down” to set time, each change with 5 minutes decrease to 0.



Red Timer for Red wavelength turn on/  
off

NIR Timer for Nir wavelength turn on/  
off



# Product Display





## WHATS INCLUDED:



LED THERAPY LIGHTS



GLASSES



COUNTRY SPECIFIC POWER CORD



DOOR BUCKLE



ADJUSTABLE HANGING KIT



DOOR HANGING KIT



DOOR HANGING KIT